




## TECH-ASSIST ALL PURPOSE SLING



**SKU:** 23411

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The All Purpose Sling is a 2-in-1 sling designed to give greater flexibility of client use. It incorporates head support to allow for a full body lift but can also be used without the head support as a G.P. style sling by folding the head section out of the way. The head section foam can also be removed for comfort if leaving behind client. Webbing attachments on sling are covered to prevent chafing. Fully washable and quick drying. S.W.L 240KG.

### SIZING

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

**SA**

**Brighton Showroom**

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

**Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

**Reynella Showroom**

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

**VIC****Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT****Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)