



TECH-ASSIST ANGLED TRANSFER BOARD

\$110.00 - \$155.00Price range: \$110.00 through

\$155.00



SKU: 16233-H

PRODUCT DESCRIPTION

A smaller transfer board angled at one end to allow better positioning around the wheel of a wheelchair for side transfers. Angled Transfer Board (P.V.C.) flexible construction for slight height variation. S.W.L 110 Kilos

FEATURES

- · A smaller transfer board
- Angled at one end to allow better positioning around the wheel of a wheelchair for side transfers
- Available in plywood which is solid, no flex and light for general transfer use
- Also available in PVC which has a flexible construction for slight height variation
- Unique low friction coating gives more slip and less grip
- Rounded edges and radius corners provide a smooth and safe transfer for both client and carer
- · Hygienic; easily cleaned