



TECH-ASSIST ANGLED TRANSFER BOARD

\$122.00 - \$172.00 Price range: \$122.00 through \$172.00



SKU: 16233-H

PRODUCT DESCRIPTION

A smaller transfer board angled at one end to allow better positioning around the wheel of a wheelchair for side transfers. Angled Transfer Board (P.V.C.) flexible construction for slight height variation.S.W.L 110 Kilos

FEATURES

- A smaller transfer board
- Angled at one end to allow better positioning around the wheel of a wheelchair for side transfers
- Available in plywood which is solid, no flex and light for general transfer use
- Also available in PVC which has a flexible construction for slight height variation
- Unique low friction coating gives more slip and less grip
- Rounded edges and radius corners provide a smooth and safe transfer for both client and carer
- Hygienic; easily cleaned