


TECH-ASSIST BANANA BOARD

\$157.00 - \$235.00 Price range: \$157.00 through \$235.00



SKU: 14339

PRODUCT DESCRIPTION

The Banana Board is used to transfer client from a vehicle seat to wheelchair, alleviating the need to lift. SWL 110KG

FEATURES

- Used to transfer people from a vehicle seat to a wheelchair
- Alleviates the need to lift
- Lightweight, smooth, flexible and resilient PVC material
- Specific shape designed to access most vehicles
- Optional grip tape on base of board
- With hand grip
- Max User Weight 110kg

SPECIFICATIONS

Max User Weight	110kg
Overall Length	1000mm
Package Quantity	1
Package Type	Single
Product Weight	1kg
Thickness	6mm