


**TECH-ASSIST FULL BODY HAMMOCK SLING**

\$437.00 - \$450.00 Price range: \$437.00 through \$450.00



SKU: 21411

PRODUCT DESCRIPTION**PURPOSE & BENEFITS**

The Full Body Hammock sling is designed to give the features of the G.P. sling but incorporating head support. This sling has a solid foam back section to make positioning behind the back easier. All foam sections used are of closed cell foam, therefore they will not absorb moisture, making the sling quick drying. S.W.L 240KG.

SIZING