



TECH-ASSIST GENERAL TRANSFER BOARD CUT OUT

\$140.00 - \$196.00 Price range: \$140.00 through \$196.00



SKU: 16133-H

PRODUCT DESCRIPTION

The general transfer board with cut out is ideal for short distance transfer from bed or chair to wheelchair, where the wheelchair has large wheel, the cut out sits around wheel and helps to alleviate sideways movement of board. Also allows for more direct transfer onto seat. SWL 110KG

FEATURES

- General transfer aid
- Suitable for a variety of short distance transfer situations, i.e. from bed to wheelchair etc.
- The cut out sits around the wheel of the wheelchair
- Alleviates sideways movement of board
- Allows for more direct transfer onto seat
- Unique low friction coating gives more slip less grip
- Rounded edges and radius corners provide a smooth and safe transfer
- Hygienic; easily cleaned
- Max User Weight 110kg