



TECH-ASSIST GENERAL TRANSFER BOARD



SKU: 15733-H

PRODUCT DESCRIPTION

The General Transfer Board as its name suggests, is used as a general transfer aid in a variety of transfer situations, ie, from bed to chair, chair to chair, chair to bed, etc. SWL 110KG

FEATURES

- General transfer aid
- Suitable for short distance transfer situations, i.e. from bed to wheelchair etc.
- Available in plywood which is a stiff lightweight transfer board, ideal in most situations
- Available in PVC which is a more flexible transfer board, ideal where there is slight height variation
- Unique low friction coating gives more slip less grip
- Rounded edges and radius corners provide a smooth and safe transfer for both client and carer
- Hygienic; easily cleaned
- Max User Weight 110kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au