




TECH-ASSIST GENERAL TRANSFER BOARD

\$124.00 - \$168.00 Price range: \$124.00 through \$168.00



SKU: 15733-H

PRODUCT DESCRIPTION

The General Transfer Board as its name suggests, is used as a general transfer aid in a variety of transfer situations, ie, from bed to chair, chair to chair, chair to bed, etc. SWL 110KG

FEATURES

- General transfer aid
- Suitable for short distance transfer situations, i.e. from bed to wheelchair etc.
- Available in plywood which is a stiff lightweight transfer board, ideal in most situations
- Available in PVC which is a more flexible transfer board, ideal where there is slight height variation
- Unique low friction coating gives more slip less grip
- Rounded edges and radius corners provide a smooth and safe transfer for both client and carer
- Hygienic; easily cleaned
- Max User Weight 110kg