




TECH-ASSIST MANUAL LIFTING SEAT



SKU: 13121-H

PRODUCT DESCRIPTION

A strong nylon or mesh seat, with four padded hand holds and a lap strap. Used to transfer client easily and securely in difficult situations. eg. wheelchair to/from car, wheelchair to/from pool or wheelchair to/from bath. Lap strap built in for extra security. Fabrics used are strong and durable, yet will not chafe the skin. Made in Australia using the highest quality materials and workmanship. Tested to International Standards, S.W.L. 150kgs* (* Please observe manual handling recommendations)

Max User Weight	150kg
Overall Depth	400mm
Overall Width	380mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au