



### **TECH-ASSIST TRANSFER BOARD WITH TWO HANDLES**

**\$151.00 - \$164.00** Price range: \$151.00 through \$164.00



**SKU:** 15433-H

## **PRODUCT DESCRIPTION**

The 6mm Transfer Board, is used as a general transfer aid in a variety of transfer situations, ie, from bed to chair, chair to chair, chair to bed, etc. A Mid Flex board, ideal where there is slight height variation. S.W.L 110 Kilos

## **FEATURES**

- Used to transfer a person from bed to chair, chair to bed etc
- Flexible
- Ideal for small gap transfers with minor differences in height
- Slide across the width of the board
- Constructed from a lightweight, smooth, resilient material
- Hygienic, can be easily wiped clean
- Rounded edges and radius corners
- Underside features grip sections ensuring secure usage