




## TECH-ASSIST TRANSLIDE MAT AND PAD



**SKU:** 15000

## PRODUCT DESCRIPTION

The Translide System is used to transfer a client who can weight-bear, from a wheelchair to the seat of a vehicle. It can also be used to help position a client, or help a client position themselves in a vehicle.

## FEATURES

- Used for transfers and positioning of clients who can weight-bear
- Non slip outer cover with single direction slip, locks in reverse
- Smooth slippery underside of pad and top cover of mat, ensures flexibility to desired position
- Able to manoeuvre in a side to side, front to back and rotationally to ensure exact positioning
- Soft foam and quilt construction provides comfort for client
- Optional rear strap on both parts to twist and secure to rear of chair
- Mat and Pad are washable
- Items sold individually or as a set

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)