




TECH-ASSIST WALKING BELT



SKU: 40851

PRODUCT DESCRIPTION

Exercise is an important part of the patients therapy, whether it be part of a rehabilitation programme, or keeping the independent client mobile. The Walk Transfer Belt was designed by Tech-Assist Products to promote walking as therapy, whilst providing comfort and security. The Walk Transfer Belt is worn around the waist, the client is able to be stabilised and guided by the carer, using the hand holds on the belt. The Walk Transfer Belt can also be used as a transfer belt for stand transfer or positional guidance. Walk/Transfer belts have optional rubber backing for extra grip to client. Walk Transfer Belts are tested to International Standards at 150 kgs.

FEATURES

- A quick action tightening cord can be pulled to provide a snug fit after the front buckle has been fastened
- Fully washable
- Is shaped on the lower edge so as not to cut into the wearers legs in a sitting position
- Optional rubber backing on nylon material belts to prevent slipping
- Various closures to choose from, Velcros, Side release buckle, Camloc buckle
- Four padded hand holds allow variation for position of carers.
- Used to assist client from sitting to standing position (or vice-versa)
- Used to stabilise and support a client whilst walking, either with or without a walking frame (depending upon clients condition)

SIZING

- Extra Small fits waist 50cm to 75cm
- Small fits waist 65 cm to 95 cm
- Medium fits waist 80 cm to 110cm
- Large fits waist 100cm to 130cm
- Extra large fits waist 125cm to 155cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au