


**TECH-ASSIST ALL PURPOSE HIGH CUT SLING****SKU: 24411****PRODUCT DESCRIPTION****PURPOSE & BENEFITS**

The High Cut Sling is essentially the same cut as a standard All Purpose sling but cut higher in the centre section to achieve more access to the waist line of client for help in undressing and toileting. Although not as accessible as the full toileting sling it certainly comes half way between the full sling and toileting. The advantage of the High Cut sling is that the client has their arms inside the sling versus the toileting sling which has support from under the arms (outside the sling). The High Cut sling incorporates a chest strap to increase support. S.W.L 240KG.

SIZING

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA**Brighton Showroom**

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au