




ALLCARE TENS MACHINE



SKU: ACPROTENS

PRODUCT DESCRIPTION

KIT OUTLINE

Kit Includes 1x TENS Machine, 4x Square Electrodes, 1x Battery, 2x Leads, 1x Carry case and 1x Instruction Booklet.

A regular TENS with adjustable pulse width, pulse rate, intensity, modes and an accurate timer of 15, 30 and 60 minutes, will stop output at pre-setting time.

NOTE: If you are unsure on how to use TENS machines correctly, we recommend seeking advice from a physician before use.

FEATURES

- Dual Channel
- Hard Carry Case
- Battery Powered
- Adjustable pulse width, pulse rate plus intensity
- Timer of 15, 30 and 60 minutes
- Will stop output at pre-set time
- Compatible with commonly available electrodes

SPECIFICATIONS

- Channel: Dual, isolated between channels
- Pulse Amplitude: 0-80mA peak to peak into 500 ohm load.
- Pulse Rate: Adjustable, from 2Hz ~ 150 Hz

- Pulse Width: Adjustable, from 30uS ~ 260 uS
- Mode: B, N, M TIMER: Selectable, 15, 30, 60 Minutes (or 15, 30, constant)
- Wave Form: Asymmetrical Bi-Phasic Square Pulse
- Voltage: 0-40V peak into 500ohm load

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au