



THERMOSMART AIRSPIRAL BREATHING TUBE F&P SLEEPSTYLE

~~\$113.00~~ Original price was: \$113.00. ~~\$83.00~~ Current price is: \$83.00. Sale Available Online Only! While Stocks Last



SKU: 900SPS120

PRODUCT DESCRIPTION

Purpose & Benefits

The **ThermoSmart AirSpiral Breathing Tube F&P SleepStyle** is a heated CPAP tubing designed exclusively for Fisher & Paykel SleepStyle Auto CPAP machines. It uses an internal heated coil to help prevent condensation ("rainout") and maintain optimal humidification levels during therapy, improving comfort and consistency throughout the night.

Features

- Designed for exclusive use with Fisher & Paykel SleepStyle Series Auto CPAP machines
- Built-in heated coil reduces condensation and supports stable humidity
- Promotes greater comfort and effective CPAP therapy

Specifications

- **Compatibility:** Fisher & Paykel SleepStyle Auto CPAP machines
- **Type:** Heated tubing