



## THRONE BARIATRIC TOILET RAIL

**\$610.00**



**SKU: RB150**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Bariatric Toilet Rail is designed to assist the larger build. This particular rail can hold weight up to 150 kg to 180 kg.

With its freestanding feet and rubber grips, the Bariatric rail relies on the floor of the bathroom/toilet, rather than being supported by the toilet itself, providing excellent stability. The legs and feet are set back to be of no hindrance to the user.

### SPECIFICATIONS

- Distance Between Handrails: Approx 600 mm
- Height Above Seat: 250 mm (with adjustable legs)
- Overall Width: 750 mm

<https://youtu.be/ti6wiFWmTwQ>