



THRONE CHILD STEP

\$295.00



SKU: STEP

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Throne Child Step is a thoughtfully designed accessory that aims to enhance bathroom safety and independence for children. It's a sturdy, user-friendly step that helps young ones access the toilet easily, promoting their autonomy and confidence in using the bathroom on their own.

FEATURES

- **Stable and Secure:** The step is designed with a non-slip surface and strong base, ensuring that it remains stable while the child uses it, reducing the risk of slips or falls.
- **Child-Friendly Height:** Tailored to bring the toilet to a comfortable height for children, making it easier for them to use it independently.
- **Durable Material:** Constructed from high-quality materials that can withstand the rigours of daily use by active children.
- **Easy to Clean:** The step's surface is easy to wipe clean, a crucial feature for maintaining hygiene in the bathroom.
- **Portable Design:** Lightweight enough to be moved and positioned by adults as needed, but sturdy enough to provide stable support for children.
- **Versatile Use:** While primarily used for toilet training and use, it can also serve as a step for reaching the sink or other areas in the bathroom.

The Throne Child Step is an excellent tool for fostering independence in young children as they learn to use the bathroom. Providing them with a safe and appropriately sized way to access the toilet and

sink, helps to build their confidence and autonomy during their formative years.

SPECIFICATIONS

- Height adjustable: 2 positions 22.5cm and 26.5cm
- Weight Capacity: 30kg
- Fits all Throne Adult and Children's Toilet Support Rails