



THRONE STANDARD 3 IN 1 TOILET RAILS



SKU: R31P0

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed for safety, comfort and benefits from the average build to the larger build persons. The 3 in 1 Rail is designed for safety and comfort and benefits a range of users from the average build to the larger build person's. This Rail enables the user to have support close to the body, where it is most needed. This model has a weight loading of 120kg.

FEATURES

- The user is able to lower safely onto the toilet and pushing up from the seat is also much easier.
- The Rail is safer than over frames as it is rigidly fixed to the toilet pan, and also does not need to be removed for other users.
- The handles can also be folded down in order to facilitate a wheelchair transfer.
- Additionally, either one or both handles can be removed to suit some requirements.

SPECIFICATIONS

Two optional settings:

- 510mm or 640mm between the handles, and fold down handles
- Height of handles above seat: 260mm approx
- The back plate internal measurements: 370mm

https://www.youtube.com/embed/hnViTA7aq6k



For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773 brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au