

THRONE STANDARD 3 IN 1 TOILET RAILS

\$499.00



SKU: R31PC



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed for safety, comfort and benefits from the average build to the larger build persons. The 3 in 1 Rail is designed for safety and comfort and benefits a range of users from the average build to the larger build person's. This Rail enables the user to have support close to the body, where it is most needed. This model has a weight loading of 120kg.

FEATURES

- The user is able to lower safely onto the toilet and pushing up from the seat is also much easier.
- The Rail is safer than over frames as it is rigidly fixed to the toilet pan, and also does not need to be removed for other users.
- The handles can also be folded down in order to facilitate a wheelchair transfer.
- Additionally, either one or both handles can be removed to suit some requirements.

SPECIFICATIONS

Two optional settings:

- 510mm or 640mm between the handles, and fold down handles
- Height of handles above seat: 260mm approx
- The back plate internal measurements: 370mm

<https://www.youtube.com/embed/hnViTA7aq6k>