


THRONE TOILET SPACER

\$130.00 - \$134.00 Price range: \$130.00 through \$134.00



SKU: SP-HEAD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Throne Toilet Spacers have been specifically designed to raise toilet seats. The spacer is supplied with longer bolts to secure the seat and spacer to the pan.

FEATURES

- The Spacers come in two sizes, 50mm and 80mm
- They can be used with or without Throne rails
- Throne spacers allow use of a normal toilet seat and they eliminate movement, which can lead to slipping off the pan.
- Spacers are easily installed, becoming an extension of the toilet pan.
- Throne Spacers are also used to raise the height of most electronic bidets.
- Spacers are easily installed, in effect becoming an extension of the toilet pan. In order to stop any lateral movement of the Spacer on the pan, use a small portion of double sided tape.
- Increases height of toilet seat

SPECIFICATIONS

- Height: 50mm or 80mm
- Product Weight: 1.5kg

<https://youtu.be/I5TaOyHOiTU>