

THRONE TOILET SPACER





SKU: SP-HEAD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Throne Toilet Spacers have been specifically designed to raise toilet seats. The spacer is supplied with longer bolts to secure the seat and spacer to the pan.

FEATURES

- The Spacers come in two sizes, 50mm and 80mm
- They can be used with or without Throne rails
- Throne spacers allow use of a normal toilet seat and they eliminate movement, which can lead to slipping off the pan.
- Spacers are easily installed, becoming an extension of the toilet pan.
- Throne Spacers are also used to raise the height of most electronic bidets.
- Spacers are easily installed, in effect becoming an extension of the toilet pan. In order to stop any lateral movement of the Spacer on the pan, use a small portion of double sided tape.
- Increases height of toilet seat



SPECIFICATIONS

Height: 50mm or 80mmProduct Weight: 1.5kg

https://youtu.be/I5TaOyHOiTU

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au