

TIDEWAVE WOUND CARE TURNING MATTRESS



SKU: 90-200-001-EU



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Tidewave Wound Care Mattress represents a significant innovation in care technology, specifically designed to assist in the regular turning of patients who are bedridden or have limited mobility. This advanced mattress system is engineered to improve patient comfort, reduce the risk of pressure ulcers, and aid caregivers in managing patient positioning with ease and efficiency.

FEATURES

- **Automated Turning:** The mattress features an automated turning mechanism that gently and periodically rotates the patient. This automated process helps in evenly distributing pressure and reducing the risk of developing pressure sores, a common concern for individuals with limited mobility.
- **Pressure Redistribution:** The mattress is constructed with high-quality materials that conform to the patient's body, providing excellent pressure redistribution across various contact points. This design minimises the risk of pressure points forming, further preventing discomfort and ulcers.
- **Quiet Operation:** Designed to operate quietly, the Tidewave Turning Mattress ensures that patients are not disturbed during sleep or rest periods when the mattress activates its turning function.
- **Ease of Use:** The mattress system is user-friendly, with intuitive controls that make it easy for caregivers to adjust settings and manage patient care effectively.
- **Safety Features:** Equipped with safety mechanisms to prevent accidents during operation, the mattress is designed with patient security in mind. These features ensure that the mattress operates smoothly and safely during each turning cycle.

- **Comfort and Support:** Beyond its turning functionality, the mattress provides a comfortable and supportive sleeping surface, catering to the overall well-being of the patient.

The Tidewave Turning Mattress is a revolutionary care solution that combines technology and comfort to significantly improve the quality of life for bedridden patients and streamline care processes for healthcare providers. By automating the turning process, it not only enhances patient care but also supports caregivers in their daily tasks, making it an invaluable addition to any care setting.

SPECIFICATIONS

- **Turning degree:** both sides, only the left or the right side
- **Turning angle:** 10, 20 or 30 degrees
- **Delay on each side:** no delay, 15 min or 30 min delay on each turning session
- **Turning speed:** low, medium or full speed*
- **Care function:** When clicking the start/pause button, the mattress deflates and starts operating automatically after 30 minutes.
- **Mattress Sizing Options:**
 - Paediatric: 700x1600 (Soft), SWL 9-40 kg
 - Adult: 900x2000 (Soft, Medium, Firm), SWL 40-150kg
 - King Single: 1070x2030 (Soft, Medium, Firm), SWL 40-150kg
 - Custom sizes available
- **Other Mattress Options**
 - **Soft** 40kg - 65kg
 - **Medium** 60kg - 115kg
 - **Firm** 110kg - 150kg

USER MANUAL

[Click here for the User Manual](#)

TIDEWAVE MATTRESS - TIMELAPSE

https://medimart.com.au/wp-content/uploads/2024/04/video_hd.mp4