




TOPFORM ASHLEY DUAL MOTOR LIFT RECLINER - MINI

\$2,415.56



SKU: 1001614

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Experience unparalleled comfort and accessibility with the Topform Ashley Lift Recliner. Equipped with dual motors, this recliner offers a unique blend of luxury and practicality, ensuring users can effortlessly transition between sitting, reclining, and standing positions. Its innovative design caters to those who value both functionality and style, making it an excellent addition to any living space.

Each chair in the Topform Recliner Ashley range is designed to suit a particular body shape and size.

FEATURES

- **Dual Motor System:** Allows for independent operation of the backrest and footrest. Users can adjust each section according to their comfort without affecting the other.
- **Lift Functionality:** Facilitates smooth transitions from a seated to a standing position, aiding those with mobility challenges.
- **Plush Upholstery:** Provides ultimate comfort, making it ideal for extended use.
- **User-friendly Controls:** Intuitive remote control for easy adjustments.
- **Sturdy Construction:** Ensures durability and longevity.
- **Integrated Side Pockets:** Convenient storage for personal items or the remote control.
- **Modern Design:** Complements a variety of home decor.

SPECIFICATIONS

Mini

- Seat Width: 48cm
- Seat Height: 49cm

- Seat Depth: 47cm
- Seat to top of back: 60cm
- Recommended user height: 153-160cm