


## TOPFORM ASHLEY SINGLE MOTOR LIFT RECLINER - MINI



**SKU:** 1001609F

### PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Topform Ashley Lift Recliner with a single motor is an elegant and functional addition to any living space, offering comfort and support for individuals who may need assistance when transitioning from a seated to a standing position. Its single-motor design allows for smooth operation and easy control of the recliner's position.

### FEATURES

- **Single Motor Operation:** Simplifies the control mechanism, allowing users to adjust the back and leg rest simultaneously with just the push of a button.
- **Lifting Function:** Gently guides users to a standing position, ideal for those who experience difficulty standing up from a seated posture.
- **Reclining Comfort:** Users can easily recline to a comfortable position for resting or watching television.
- **Quality Upholstery:** The recliner comes in a variety of fabric or leather choices to fit any décor while providing a comfortable and durable seating surface.
- **Battery Backup:** Ensures functionality even in the event of a power outage, providing peace of mind and increased safety.
- **Sturdy Construction:** Built to accommodate individuals of various sizes, offering reliable stability and long-lasting durability.
- **Convenient Size:** Designed to fit well in both small and large spaces, providing the benefits of a lift chair without overwhelming the room.

The Topform Ashley Lift Recliner with a single motor provides an optimal blend of comfort, functionality, and style, designed to meet the needs of individuals requiring mobility assistance. With

its ease of use and reliable performance, this lift recliner is a popular choice for those looking to improve their quality of life at home.

## SPECIFICATIONS

### Mini

- Seat Width: 48cm
- Seat Height: 49cm
- Seat Depth: 47cm
- Seat to top of back: 60cm
- Recommended user height: 153-160cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)