



TOPRO NEURO SAFETY BRAKE WALKER WITH BACKREST



SKU: 814746-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Topro Neuro Safety Brake Walker is designed for people with Parkinson's disease or other neurological conditions. It features a reverse braking system that automatically locks movement unless the braking levers are compressed. This can help to prevent falling and ensure safety when sitting down and standing up.

The rollator also promotes good posture by encouraging users to walk inside the frame, rather than behind it. This can help to improve balance and reduce the risk of falls. Additionally, the rollator folds compactly and can be transported in cars, planes, trains, and buses.

FEATURES

- The reverse braking system automatically locks movement unless the braking levers compressed
- Promotes good posture by walking inside the frame, not behind it
- Backrest included
- Available in Small & Medium
- Folds compact and can travel in cars, planes, trains & buses

[Check out the Drag Brake attachment here](#)

Why you need it

If you have Parkinson's disease or other neurological conditions, the Topro Neuro Safety Brake Walker can help you stay safe and mobile. The reverse braking system can help to prevent falls, the backrest can provide support, and the compact folding design makes it easy to transport.

Order your TOPRO Troja Neuro rollator today!

SPECIFICATIONS

	Medium	Small
Colour	Dark Grey	Silver
Recommended user height	150-200 cm	135-170 cm
Maximum weight of user	150 kg	125 kg
Capacity shopping bag	10 kg / 16 l	10 kg / 16 l
Handle height	78-100 cm	67-86 cm
Maximum length	65 cm	65 cm
Maximum width	60 cm	60 cm
Width between handles	47 cm	47 cm
Handles – Ergonomic, different holding positions (TOPRO Ergo Grip)	Handles – Ergonomic, different holding positions (TOPRO Ergo Grip)	Handles – Ergonomic, different holding positions (TOPRO Ergo Grip)
Seat height	62 cm	54 cm
Diameter turning circle	84 cm	84 cm
Diameter wheels	20 cm	20 cm
Width wheels	3.6 cm	3.6 cm
Grip distance driving brake	6.5 cm	6.5 cm
Height when folded	80 cm	74 cm
Width when folded	24 cm	24 cm
Length when folded	65 cm	65 cm
Product weight STANDARD model	7.8 kg	7.4 kg
Art. No. STANDARD model	814746	814747

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au