



TOPRO OLYMPOS ATR OFF-ROAD ROLLATOR WITH BACKREST



SKU: 814300

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Topro Olympos ATR Off-Road Rollator with Backrest is a premium mobility aid designed to handle rough terrains while providing supreme comfort and stability. This rollator is ideal for users who love outdoor adventures but require additional support.

FEATURES

- **Rugged Design for Off-Road Use:** The Olympos ATR is built to conquer challenging terrains. Its large, sturdy wheels ensure smooth movement across uneven surfaces, making it perfect for outdoor excursions.
- **Ergonomic Handles:** The rollator features adjustable, ergonomic handles that provide a comfortable grip, reducing strain on the hands and wrists.
- **Built-in Seat with Backrest:** For moments when you need to take a break, the Olympos ATR offers a convenient and comfortable seat, complete with a supportive backrest.
- **Secure Braking System:** Equipped with an efficient and reliable braking system, this rollator ensures your safety, giving you complete control over your movements.
- **Easy to Fold and Transport:** The rollator is designed for easy folding, making it compact and simple to transport or store when not in use.
- **Sturdy Frame:** The robust frame of the Olympos ATR guarantees long-lasting durability, providing reliable support for users.
- **Spacious Storage:** With its built-in storage solutions, you can easily carry your essentials with you on your outdoor adventures.

SPECIFICATIONS

- **Dimensions:**
 - Folded: 77 x 27 x 87.5 cm
 - Unfolded: 77 x 68 x 97 cm
- **Weight:** 11kg
- **Weight capacity:** 136kg
- **Seat height:** 60 cm
- **Handle height:** 78.5-94.5 cm (31-37 inches)
- **Wheels:**
 - Front: 20 cm (8 inches)
 - Rear: 15 cm (6 inches)

BROCHURE

[Topro Rollators Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au