




## TOPRO TROJA 5G ROLLATOR



**SKU:** 814600-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Topro Troja 5G Rollator represents a pinnacle of design and functionality in the world of mobility aids. This state-of-the-art walking support is crafted for individuals who desire a combination of stability, style, and innovation as they navigate various terrains.

### FEATURES

- **Advanced Braking System:** The Troja 5G features an intuitive and reliable braking system, providing users with confidence and security when moving or coming to a stop.
- **Lightweight Construction:** Manufactured from high-quality materials, the rollator is light yet sturdy, allowing for easy transportation and manoeuvrability without sacrificing durability.
- **Ergonomic Handles:** Ergonomically shaped handles offer a comfortable grip, aiding in user comfort and reducing hand strain during prolonged use.
- **Flexible Seating Option:** It comes equipped with a built-in seat, giving users a convenient place to rest when needed.
- **Adjustable Height:** The handles are fully adjustable, ensuring users of various heights can find a comfortable and ergonomic fit.
- **Foldable Design:** The Troja 5G can be easily folded, making it compact for storage or fitting into a car boot when travelling.
- **Accessory Compatibility:** It is designed to be compatible with a range of accessories, such as baskets, trays, and cane holders, allowing users to customise the rollator to their personal

needs.

The Topro Troja 5G Rollator is not only a mobility tool but also a statement of independence and freedom, helping users maintain an active and fulfilling lifestyle. With its innovative features and ergonomic design, it has become a popular choice for those seeking a high-end mobility solution.

## SPECIFICATIONS

- Weight: 22.6 lbs (10.2 kg)
- Folded dimensions: 28.3" x 19.7" x 15.7" (72 cm x 50 cm x 40 cm)
- Seat height: 22.8" - 31.5" (58 - 80 cm)
- Handle height: 35.4" - 42.9" (90 - 109 cm)
- Wheel size: 8" (20 cm)

## BROCHURE

[Topro Rollators Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)