



# **TOPRO TROJA 5G ROLLATOR**

\$658.00 - \$723.00Price range: \$658.00 through \$723.00



**SKU:** 814600-H

#### PRODUCT DESCRIPTION

### **PURPOSE & BENEFITS**

The Topro Troja 5G Rollator represents a pinnacle of design and functionality in the world of mobility aids. This state-of-the-art walking support is crafted for individuals who desire a combination of stability, style, and innovation as they navigate various terrains.

## **FEATURES**

- **Advanced Braking System:** The Troja 5G features an intuitive and reliable braking system, providing users with confidence and security when moving or coming to a stop.
- **Lightweight Construction:** Manufactured from high-quality materials, the rollator is light yet sturdy, allowing for easy transportation and manoeuvrability without sacrificing durability.
- **Ergonomic Handles:** Ergonomically shaped handles offer a comfortable grip, aiding in user comfort and reducing hand strain during prolonged use.
- Flexible Seating Option: It comes equipped with a built-in seat, giving users a convenient place to rest when needed.
- **Adjustable Height:** The handles are fully adjustable, ensuring users of various heights can find a comfortable and ergonomic fit.
- **Foldable Design:** The Troja 5G can be easily folded, making it compact for storage or fitting into a car boot when travelling.
- **Accessory Compatibility:** It is designed to be compatible with a range of accessories, such as baskets, trays, and cane holders, allowing users to customise the rollator to their personal



needs.

The Topro Troja 5G Rollator is not only a mobility tool but also a statement of independence and freedom, helping users maintain an active and fulfilling lifestyle. With its innovative features and ergonomic design, it has become a popular choice for those seeking a high-end mobility solution.

### **SPECIFICATIONS**

• Weight: 22.6 lbs (10.2 kg)

• Folded dimensions: 28.3" x 19.7" x 15.7" (72 cm x 50 cm x 40 cm)

• Seat height: 22.8" - 31.5" (58 - 80 cm)

• Handle height: 35.4" - 42.9" (90 - 109 cm)

• Wheel size: 8" (20 cm)

#### **BROCHURE**

**Topro Rollators Brochure**