



## TOPRO TROJA WALKER WITH FOREARM GUTTER



**SKU:** 815430

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Topro Troja Walker with Forearm Gutter is a sturdy and lightweight walker designed to provide support and stability for people with mobility challenges. The forearm gutter design allows you to walk inside the frame, which encourages good posture and reduces fatigue. The walker also features a comfortable padded seat, adjustable handlebars, and brakes on both wheels for added safety.

### FEATURES

- Lightweight and portable design
- Adjustable handlebar height for a comfortable fit
- Forearm gutter design for good posture and reduced fatigue
- Padded seat for resting
- Brakes on both wheels for added safety
- Folds easily for storage and transportation

The Topro Troja Walker with Forearm Gutter is a great option for people who need a stable and supportive walker for indoor use. It is available in three sizes to fit a variety of users.

### SPECIFICATIONS

#### Size M (medium)

- **Recommended user height:** 150-185 cm
- **Max. user weight:** 135 kg
- **Max. capacity basket:** 5 kg / 20 l
- **Height of forearm support points:** 98-115 cm

- **Handles:** Round, Ø 3.5 cm
- **Maximum length:** 65 cm
- **Maximum width:** 63 cm
- **Seat height:** 63 cm
- **Seat width:** 45,5 cm
- **Turning circle:** 84 cm
- **Diameter wheels:** 20 cm
- **Width wheels:** 3.6 cm
- **Grip distance driving brake:** 7.4 cm
- **Height when folded:** 102 cm
- **Width when folded:** 33 cm
- **Length when folded:** 65 cm
- **Total dimensions:** (l × w × h) 65 × 33 × 102 cm
- **Product weight:** 8.9 kg

## BROCHURE

[Topro Rollators Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)