



TOPRO TROJA WALKER WITH FOREARM GUTTER

\$1,427.00



SKU: 815430

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Topro Troja Walker with Forearm Gutter is a sturdy and lightweight walker designed to provide support and stability for people with mobility challenges. The forearm gutter design allows you to walk inside the frame, which encourages good posture and reduces fatigue. The walker also features a comfortable padded seat, adjustable handlebars, and brakes on both wheels for added safety.

FEATURES

- Lightweight and portable design
- Adjustable handlebar height for a comfortable fit
- Forearm gutter design for good posture and reduced fatigue
- Padded seat for resting
- Brakes on both wheels for added safety
- Folds easily for storage and transportation

The Topro Troja Walker with Forearm Gutter is a great option for people who need a stable and supportive walker for indoor use. It is available in three sizes to fit a variety of users.

SPECIFICATIONS

Size M (medium)

- **Recommended user height:** 150-185 cm
- **Max. user weight:** 135 kg
- **Max. capacity basket:** 5 kg / 20 l
- **Height of forearm support points:** 98-115 cm

- **Handles:** Round, Ø 3.5 cm
- **Maximum length:** 65 cm
- **Maximum width:** 63 cm
- **Seat height:** 63 cm
- **Seat width:** 45,5 cm
- **Turning circle:** 84 cm
- **Diameter wheels:** 20 cm
- **Width wheels:** 3.6 cm
- **Grip distance driving brake:** 7.4 cm
- **Height when folded:** 102 cm
- **Width when folded:** 33 cm
- **Length when folded:** 65 cm
- **Total dimensions:** (l × w × h) 65 × 33 × 102 cm
- **Product weight:** 8.9 kg

BROCHURE

[Topro Rollators Brochure](#)