



TRIVIDA 24" WHEELCHAIR WHEEL SET



SKU: CC20-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Trivida Detachable Wheelchair Wheel is a groundbreaking invention that aims to improve the daily mobility of wheelchair users. Understanding the challenges faced by those using wheelchairs, Trivida's solution seeks to simplify transfers and promote a barrier-free life.

FEATURES

- **Easier Transfers:** These wheels are strategically designed to facilitate seamless transfers. Whether you're shifting from your wheelchair to the bed, car, or any other platform, Trivida ensures that the process is hassle-free.
- **Promote Day-to-Day Mobility:** Wheelchair users often face challenges in their day-to-day lives when it comes to transfers. The Trivida wheel aims to reduce the energy required for such transfers, making everyday tasks simpler and more manageable.
- **Empowerment & Independence:** For many wheelchair users, relying on assistance for transfers can be a challenging experience, both mentally and physically. The Trivida wheel encourages self-reliance, enabling users to handle transfers on their own, and fostering a sense of independence.
- **User-Friendly Design:** The steps to use the detachable wheel are straightforward:
 - **Position:** Ensure the wheelchair is in the right place for a transfer.
 - **Brake:** Apply the brakes to ensure the wheelchair remains stationary.
 - **Detach:** The top segment of the wheel can easily be removed.
 - **Transfer:** With the wheel segment removed, users can safely and easily transfer.
- **Beneficial for Carers:** It's not just the users; even their family members or caregivers find the task of transferring a user exhausting. Trivida's design ensures that transfers are simpler,

requiring less physical exertion.

- **A Solution for All Ages:** As wheelchair users age, they may face a reduction in strength and mobility. The Trivida wheel accounts for these challenges, ensuring that users of all ages can benefit from its design.

The Trivida Detachable Wheelchair Wheel is more than just a product; it's a step towards a more inclusive and accessible world. By addressing one of the primary challenges faced by wheelchair users, paves the way for a life with fewer barriers and greater mobility.

SPECIFICATIONS

- Available in: Plastic PT20, Carbon CC20
- Wheel diameter: 24"
- Maximum permissible camber: 3°
- Load-bearing capacity: Wheelchair 150 kg
- Axle Quick-change axle
- 2 bearings for enhanced stability
- Weight of wheel: 3900 g Plastic, 3300 g Carbon
- Design: Spoke wheel, matt surface, black High-gloss | carbon disc, black

BROCHURE

[Trivida Wheel Set Brochure](#)

[Trivida Wheel Set Information Sheet](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au