




## **U-STEP WALKER 2 WITH OPTIONAL LASER AND SOUND CUE**



**SKU: USTEP2**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The U-Step Walker 2 has been created to increase your independence. The secret is in our patented U-shaped base. Its ultra-stable foundation braces you in every direction.

It is not like pushing a typical walker. Instead, the U-Step 2 surrounds you and moves with you. You will feel as stable as you would feel holding onto another person's arm.

The innovative braking system is easy to use and puts you in complete control. The U-Step 2 will not roll unless you are ready to walk. When you lightly squeeze either brake release lever, the unit will roll with you. Once you release the lever, the unit will stop immediately. This feature is particularly helpful when standing up from a chair because the unit will not roll away from you.

### **KEY FEATURES**

- Reversed braking system
- Small turning circle for tight spaces
- Rolling resistance control
- Spring-loaded front wheels
- Folds for easy transporting
- Laser & sound cue module attachments are available

### **LASER AND SOUND CUE**

This optional feature is primarily used by those with Parkinson's but can be used by anyone with an irregular gait pattern.

This Cueing Module can help get you started, normalise your walking and increase your stride. Press the red button on the module to project a bright red laser line on the floor to guide your step. You can also turn on the sound feature to set a beat pattern for walking speed.

## PACKAGES

- To purchase both the U-Step Walker 2 and Laser + Sound Cue option - Please add both the 'U-Step 2 Walker 2 Only' and 'Laser & Sound Cue Only' to your cart.
  - A \$100 service charge may apply for the manufacturer assembly of fitting the laser to the U-Step walker\*
- To purchase just the U-Step Walker 2 - Please add the 'U-Step 2 Walker Only' to your cart.

## MORE INFO

[embed]https://youtu.be/oEimu30rm4s[/embed]

Buy this U Step Walker online at Medimart as we have NDIS Payment Method plus AfterPay and ZipPay.

## SPECIFICATIONS

Overall width	58cm
Seat height	56cm
Length	63.5cm
Turning circle	74cm
Seat width	51cm
Turning circle	74cm
Weight with seat and basket	10kg
Weight capacity	170kg
Folded dimensions	107cm x 58cm x 25cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North Brighton  
 SA 5048  
 08 7079 8773  
[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave,  
 Elizabeth Vale SA 5112  
 08 7079 8777  
[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella  
 SA 5161  
 08 7079 8771  
[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)