



U-STEP WALKER 2 WITH OPTIONAL LASER AND SOUND CUE

\$685.00 - \$2,075.00 Price range: \$685.00 through \$2,075.00



SKU: USTEP2

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The U-Step Walker 2 has been created to increase your independence. The secret is in our patented U-shaped base. Its ultra-stable foundation braces you in every direction.

It is not like pushing a typical walker. Instead, the U-Step 2 surrounds you and moves with you. You will feel as stable as you would feel holding onto another person's arm.

The innovative braking system is easy to use and puts you in complete control. The U-Step 2 will not roll unless you are ready to walk. When you lightly squeeze either brake release lever, the unit will roll with you. Once you release the lever, the unit will stop immediately. This feature is particularly helpful when standing up from a chair because the unit will not roll away from you.

KEY FEATURES

- Reversed braking system
- Small turning circle for tight spaces
- Rolling resistance control
- Spring-loaded front wheels
- Folds for easy transporting
- Laser & sound cue module attachments are available

LASER AND SOUND CUE

This optional feature is primarily used by those with Parkinson's but can be used by anyone with an irregular gait pattern.

This Cueing Module can help get you started, normalise your walking and increase your stride. Press the red button on the module to project a bright red laser line on the floor to guide your step. You can also turn on the sound feature to set a beat pattern for walking speed.

PACKAGES

- To purchase both the U-Step Walker 2 and Laser + Sound Cue option - Please add both the 'U-Step 2 Walker 2 Only' and 'Laser & Sound Cue Only' to your cart.
 - A \$100 service charge may apply for the manufacturer assembly of fitting the laser to the U-Step walker*
- To purchase just the U-Step Walker 2 - Please add the 'U-Step 2 Walker Only' to your cart.

MORE INFO

[embed]<https://youtu.be/oEimu30rm4s>[/embed]

Buy this U Step Walker online at Medimart as we have NDIS Payment Method plus AfterPay and ZipPay.

SPECIFICATIONS

| | |
|-----------------------------|---------------------|
| Overall width | 58cm |
| Seat height | 56cm |
| Length | 63.5cm |
| Turning circle | 74cm |
| Seat width | 51cm |
| Turning circle | 74cm |
| Weight with seat and basket | 10kg |
| Weight capacity | 170kg |
| Folded dimensions | 107cm x 58cm x 25cm |