



## URIMED CATHETER 2 WAY 14 CH UNISEX 10ML SILICONE 40CM

**\$18.00** Original price was: \$18.00.<u>\$16.00</u>Current price is: \$16.00.



**SKU:** 10002156

## **PRODUCT DESCRIPTION**

## **FEATURES AND BENEFITS**

Urimed catheters are 100% silicone material which is very soft and hypoallergenic. Its 40cm length is suitable for both male and female.

It has a radiopaque tip and shaft with a short bullet tip to avoid irritation of the bladder. Latex and PVC free. Available in CH12 - CH24.