



URIBAG PORTABLE URINAL FOR MALES



SKU: SYMURIBAG-M

PRODUCT DESCRIPTION

PORTABLE URINAL PURPOSE AND BENEFITS

The UriBag Male Travel Portable Urinal is a convenient solution for wheelchair users, men who are incontinent and those needing to urinate where there are no convenient facilities. The compact storage tube, efficient shape and Latex construction ensure the Uribag male urinal is lightweight and durable. The snap-lock seal does not leak, allowing it to be discreetly used, sealed, then stored for emptying and washing later. Available whenever needed, the UriBag can be carried in a pocket, bag or glove box and requires only minimal privacy for use. Suitable not only for those with a disability but also as a practical solution for travelling or camping.

KEY FEATURES

- Compact modern design, stores privately in pockets or bags when not in use.
- In-built storage tube to easily tuck away the latex bladder.
- Odour proof for maximum discretion.
- Leak proof lid, seals Uribag to prevent spills. Must be attached before setting bag down.
- Suitable for use whilst sitting, standing or lying down.

SPECIFICATIONS

Folded height 8cm

Folded diameter 4.5cm

Maximum fluid capacity 1200ml

Recommended replacement frequency 6-12 months (with normal use).

Available delivered Australia wide.

Also, check [spill proof male urinal](#).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au