



URIBAG PORTABLE URINAL FOR MALES

\$73.00



SKU: SYMURIBAG-M

PRODUCT DESCRIPTION

PORTABLE URINAL PURPOSE AND BENEFITS

The UriBag Male Travel Portable Urinal is a convenient solution for wheelchair users, men who are incontinent and those needing to urinate where there are no convenient facilities. The compact storage tube, efficient shape and Latex construction ensure the Uribag male urinal is lightweight and durable. The snap-lock seal does not leak, allowing it to be discreetly used, sealed, then stored for emptying and washing later. Available whenever needed, the UriBag can be carried in a pocket, bag or glove box and requires only minimal privacy for use. Suitable not only for those with a disability but also as a practical solution for travelling or camping.

KEY FEATURES

- Compact modern design, stores privately in pockets or bags when not in use.
- In-built storage tube to easily tuck away the latex bladder.
- Odour proof for maximum discretion.
- Leak proof lid, seals Uribag to prevent spills. Must be attached before setting bag down.
- Suitable for use whilst sitting, standing or lying down.

SPECIFICATIONS

Folded height 8cm

Folded diameter 4.5cm

Maximum fluid capacity 1200ml

Recommended replacement frequency 6-12 months (with normal use).

Available delivered Australia wide.

Also, check [spill proof male urinal](#).