



VARILITE 4 POINT HIP BELT



SKU: 93225-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The low-slip, laminated closed-cell foam used to make the pads is edge-bonded for comfort and longevity. The user can choose between a static and dynamic harness with the use of a dynamic positioning strap. Dynamic movement is advantageous for clients with spasticity or when doing functional tasks like reaching. The shoulder/back level is where the adjustment is made. When positioning the user from behind, carers should utilise a top-pull chest harness. When adjusting the harness, the carer gains more leverage and keeps appropriate body mechanics. There is also a bottom pull attachment available.

FEATURES

- As for 2 points with 2 additional adjustable front straps to assist keep hip belts in place and accommodate a difficult-to-position pelvis
- 40 mm webbing with options for different buckles and padding sizes
- The design allows padding to be positioned to fit each individual and secured in place with tri-glides
- Contoured shape laminated foam and binded edges
- To accommodate growth, replace padding with a larger size
- Available sizes XS, S, M, L

SPECIFICATIONS

Size	Width	Length	Height
X Small	4.5cm	117cm	1cm

Small	5.5cm	130cm	1cm
Medium	5.5cm	142cm	1cm
Large	7cm	155cm	1.5cm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au