



VARILITE 4 POINT HIP BELT

\$249.00 Online Only



SKU: 93225-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The low-slip, laminated closed-cell foam used to make the pads is edge-bonded for comfort and longevity. The user can choose between a static and dynamic harness with the use of a dynamic positioning strap. Dynamic movement is advantageous for clients with spasticity or when doing functional tasks like reaching. The shoulder/back level is where the adjustment is made. When positioning the user from behind, carers should utilise a top-pull chest harness. When adjusting the harness, the carer gains more leverage and keeps appropriate body mechanics. There is also a bottom pull attachment available.

FEATURES

- As for 2 points with 2 additional adjustable front straps to assist keep hip belts in place and accommodate a difficult-to-position pelvis
- 40 mm webbing with options for different buckles and padding sizes
- The design allows padding to be positioned to fit each individual and secured in place with tri-glides
- Contoured shape laminated foam and binded edges
- To accommodate growth, replace padding with a larger size
- Available sizes XS, S, M, L

SPECIFICATIONS

Size	Width	Length	Height
X Small	4.5cm	117cm	1cm

Small	5.5cm	130cm	1cm
Medium	5.5cm	142cm	1cm
Large	7cm	155cm	1.5cm