

VARILITE EVOLUTION WHEELCHAIR CUSHION





SKU: 76625-H

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS OF THE VARILITE EVOLUTION WHEELCHAIR CUSHION

The Varilite Evolution Wheelchair Cushion is designed to provide unmatched comfort and pressure management for wheelchair users. This state-of-the-art cushion combines air-foam floatation technology with adjustable features to cater to individual needs, enhancing overall seating posture and reducing the risk of pressure sores. Ideal for both active and sedentary users, the Varilite Evolution Wheelchair Cushion promotes stability, supports pelvic positioning, and improves user comfort during long periods of sitting. It's a perfect solution for those seeking to enhance their wheelchair's comfort without compromising on mobility or independence.

FEATURES

- **Air-Foam Floatation Technology**: Utilises a unique combination of air and foam to distribute weight evenly, minimizing pressure points and enhancing comfort for extended periods.
- **Adjustable Support**: Features an easily adjustable valve system that allows users to customize the cushion's firmness according to their comfort and support needs.
- **Moisture-Resistant Cover**: Comes with a breathable, moisture-resistant cover that promotes air circulation, keeping the user cool and dry throughout the day.
- **Lightweight Design**: The cushion's lightweight construction ensures it does not add unnecessary weight to the wheelchair, making it easier to maneuver and transport.
- **Durable Construction**: Built to withstand daily use, the cushion's materials are chosen for their durability and longevity, ensuring it maintains its shape and support over time.
- **Easy to Clean**: The cushion and its cover are designed for easy cleaning, a crucial feature for maintaining hygiene and comfort.



SPECIFICATIONS

• Cover Material: Available in Incontinence or Mesh

Width x depth (Inches)	Width x depth (cm)	Weight capacity (kg)
12x14	30x35	180
14×14	35x35	200
14x16	35x40	200
14x18	35x45	200
14x20	35×50	200
15×15	38x38	225
15×17	38x43	225
16×14	40x35	225
16×16	40×40	225
16×18	40×45	225
16×20	40×50	225
17×17	43x43	250
18×16	45×40	300
18×18	45×45	300
18x20	45x50	300
20x16	50×40	350
20x18	50x45	350
20x20	50x50	350
22x18	55x45	350
22x20	55x50	350
24x18	60x45	350
24x20	60x50	350

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048
08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777 elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771 reynella@medimart.com.au



VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au