



# VARILITE ICON PAL SWING-AWAY LATERAL SUPPORT



SKU: 89130C-H

#### PRODUCT DESCRIPTION

#### **PURPOSE AND BENEFITS**

Varilite's Positive Action Lateral, PAL™ Swing-away laterals are designed to use with the Icon Mid and Icon Tall backs, providing lateral positioning and support wherever needed. A PAL lateral system is made up of a lateral pad, shell bracket, pad bracket, and push-button hinge. A comfortable and supportive fit that improves the user's posture and overall seated experience is guaranteed by the PAL lateral system's extensive adjustment and customization options.

#### **FEATURES**

- The pads come in three sizes and two different styles: contoured or linear.
- They are constructed of closed-cell foam with a steel insert that enables the adjustability of the horizontal pad.
- Easily removable and washable cover shields the pad foam.
- An easily activated push button unlocks the PAL swing-away hinge mechanism.
- Available in Contoured or Linear design

#### **SPECIFICATIONS**

PAL Swing away lateral pad bracket Anterior Pad Adjustability

Fixed lateral pad bracket

Long 5.7cm (2.2")

For further information go to **medimart.com.au** or give us a call: **1300 633 063** 



#### SA

#### **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

#### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

#### **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

### **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

## NT

#### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

#### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au