



## VARILITE ICON PAL SWING-AWAY LATERAL SUPPORT

**\$401.00**



**SKU:** 89130C-H

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

Varilite's Positive Action Lateral, PAL™ Swing-away laterals are designed to use with the Icon Mid and Icon Tall backs, providing lateral positioning and support wherever needed. A PAL lateral system is made up of a lateral pad, shell bracket, pad bracket, and push-button hinge. A comfortable and supportive fit that improves the user's posture and overall seated experience is guaranteed by the PAL lateral system's extensive adjustment and customization options.

### FEATURES

- The pads come in three sizes and two different styles: contoured or linear.
- They are constructed of closed-cell foam with a steel insert that enables the adjustability of the horizontal pad.
- Easily removable and washable cover shields the pad foam.
- An easily activated push button unlocks the PAL swing-away hinge mechanism.
- Available in Contoured or Linear design

### SPECIFICATIONS

PAL Swing away lateral pad bracket  
 Fixed lateral pad bracket  
 Long

Anterior Pad Adjustability  
  
 5.7cm (2.2")