




## VECTOR PULSE - BLUETOOTH TENS & EMS MACHINE

**\$189.00**



**SKU: VPULSE-2**



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Vector PULSE is our leading TENS & EMS Muscle Stimulator device. The Vector PULSE is able to connect to your phone via Bluetooth and our interactive app allowing you to customise your treatment to your exact needs.

### MUSCLE STIMULATION, MINUS THE GUESSWORK

The accuracy of pad placement, intensity level, and session time of any muscle stimulation is critically important for safety and efficacy. That's why the Vector PULSE allows for every setting to be customised, all to help find the precise combination for you.

### BEHIND THE TECHNOLOGY: WHAT IS TENS AND NMES?

Transcutaneous electrical nerve stimulation (TENS) and neuromuscular electrical stimulation (NMES) send electronic pulses to stimulate your motor nerves. This is done by creating muscle contraction and stimulating sensory nerves to temporarily relieve pain.

### THE POWERFUL BENEFITS OF MUSCLE STIMULATION

**Improve Performance** - PULSE, through its precise electrical stimulation, effectively engages both slow and fast twitch muscle fibres, inducing contractions that closely replicate the effects of traditional training methods, ultimately leading to enhanced strength and improved overall performance.

**Natural Pain Relief** - TENS, with its targeted stimulation of sensory nerves, not only provides temporary pain relief but also triggers the release of endorphins, the body's natural painkillers,

offering a holistic approach to managing discomfort. Additionally, NMES, by improving circulation, aids in the elimination of metabolic waste, leading to reduced muscle soreness and contributing to a more effective pain relief experience.

**Speeds Up Recovery** - Smart recovery settings combine TENS and NMES to work together to reduce soreness and restore the performance of the muscle. NMES contracts the muscle, increasing the delivery of oxygenated, nutrient-rich blood to kickstart the recovery process.

**THE SMARTER THE TECHNOLOGY, THE SAFER AND MORE EFFECTIVE THE TREATMENT.**

## **WHAT'S INCLUDED**

- Vector PULSE
- 1x Large Electrode Pad
- 2x Small Electrode Pad
- 1x Magnetic Connector Wire
- Vector PULSE Charger USB Cable

**VECTOR PULSE REQUIRES A SMART PHONE TO WORK.**