

## VICAIR ACTIVE O2 WHEELCHAIR CUSHION

**\$954.00 - \$1,093.00** Price range: \$954.00 through \$1,093.00



**SKU:** OCT



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Vicair Active O2 is designed for the active wheelchair user. It offers a unique combination of ergonomically shaped foam and Vicair SmartCells. The four rear SmartCell compartments provide a high level of comfort and skin protection and offload the sitting bones and tailbone area. The foam front makes transfers easy and feels steady and improves wheelchair handling. An easy-to-clean fluid-resistant liner protects the foam inlay. The rest of the cushion is 100% machine washable and breathable for optimal hygiene, temperature, and moisture regulation. The Vicair Active O2 has an easily accessible storage pouch for precious belongings (keys, phone, wallet). In short: everything you need to maintain your active lifestyle!

The wheelchair cushion Vicair Active O2 helps prevent pressure ulcers. It is also a high-quality addition to a pressure ulcer treatment plan (stage I-IV).

### FEATURES

- Optimal temperature and moisture regulation
- Perfect for users who transfer frequently
- Excellent lateral and frontal stability
- Features a nifty storage pouch to have personal items in reach
- Depth: 9cm
- SWL: 200kg

### SPECIFICATIONS

Available height	6 and 9 cm
------------------	------------

Average weight	1250 grams
Max. user weight	200 kg
Ready for use	Yes
Standard sizes	145 ( <a href="#">Size Chart</a> )
Machine washable	Yes. Max. 60°C (without foam)
Drying time	Air dry 5 to 6 hours
Latex free	Yes
Warranty	24 months