


VICAIR ACTIVE O2 WHEELCHAIR CUSHION



SKU: OCT

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Vicair Active O2 is designed for the active wheelchair user. It offers a unique combination of ergonomically shaped foam and Vicair SmartCells. The four rear SmartCell compartments provide a high level of comfort and skin protection and offload the sitting bones and tailbone area. The foam front makes transfers easy and feels steady and improves wheelchair handling. An easy-to-clean fluid-resistant liner protects the foam inlay. The rest of the cushion is 100% machine washable and breathable for optimal hygiene, temperature, and moisture regulation. The Vicair Active O2 has an easily accessible storage pouch for precious belongings (keys, phone, wallet). In short: everything you need to maintain your active lifestyle!

The wheelchair cushion Vicair Active O2 helps prevent pressure ulcers. It is also a high-quality addition to a pressure ulcer treatment plan (stage I-IV).

FEATURES

- Optimal temperature and moisture regulation
- Perfect for users who transfer frequently
- Excellent lateral and frontal stability
- Features a nifty storage pouch to have personal items in reach
- Depth: 9cm
- SWL: 200kg

SPECIFICATIONS

Available height

6 and 9 cm

Average weight	1250 grams
Max. user weight	200 kg
Ready for use	Yes
Standard sizes	145 (Size Chart)
Machine washable	Yes. Max. 60°C (without foam)
Drying time	Air dry 5 to 6 hours
Latex free	Yes
Warranty	24 months

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au