

VICAIR ADJUSTER O2 WHEELCHAIR CUSHION



SKU: OAD

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Vicair Adjuster O2 is extremely suitable for use with asymmetric body contours and amputations. It automatically manages small obliquities and is individually adjustable to support more extreme asymmetry. When you're looking for a high comfort level, optimal skin protection and stability: the Adjuster O2 is always a safe choice. The design allows total immersion to equally redistribute pressure and prevents sliding forward. 100% machine washable and breathable, for optimal hygiene and perfect temperature and moisture regulation. As a result, the Vicair Adjuster O2 offers you optimal protection against skin breakdown.

The Vicair Adjuster O2 is a wheelchair cushion that helps prevent pressure ulcers. It is also a high-quality addition to a pressure ulcer treatment plan (stage I-IV).

FEATURES

- Optimal temperature and moisture regulation
- Fully machine washable
- Individual solution for asymmetry
- Minimises sliding and shear forces
- Optimal adjustment to the body's contours for equal pressure distribution
- SWL: 250kg

SPECIFICATIONS

Available height	6 and 10 cm
Average weight	725 grams – depending on size

Max. user weight	250 kg
Ready for use	Yes
Standard sizes	145 (Size Chart)
Machine washable	Yes. Max. 60°C
Drying time	Air dry 5 to 6 hours
Latex free	Yes
Warranty	24 Months

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au