

VICAIR ALLROUNDER O2 CUSHION



SKU: OAL/



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Active as you are, not all of your activities require or allow the use of a wheelchair. In these situations the activity cushion Vicair AllRounder O2 keeps you comfortable and helps to protect your skin. As it is easily buckled onto your body, it follows you around wherever you go.

Feel protected and secure during your sports- and outdoor activities, when you sit in a normal chair or on any hard surface, or when you play on the ground. Afterwards you can simply put the entire AllRounder O2 in the washing machine to clean it up and make it ready for your next adventure together. Experience the freedom!

FEATURES

- Optimal temperature and moisture regulation
- Fully machine washable
- Stays with you during activities outside of the wheelchair
- Skin protection
- Extremely light to support optimal mobility

CHOOSE THE RIGHT SIZE

The activity cushion Vicair AllRounder O2 is available in three sizes: small, medium and large and has an average weight of only 800 grams.

Measure the width of your hips to check which size matches the size of your body. When seated, hold a tape measure horizontally above the widest point of your upper legs, closely against the front of

your hips. Do not follow your leg contours, keep the tape measure straight. Measure from A to B (or B to A). Choose the larger size if the measured size is in between two sizes.

Small < 38 cm / < 15 inch

Medium 38 – 42,5 cm / 15- 16¾ inch

Large 43 – 48 cm / 17 – 19 inch

SPECIFICATIONS

Height	n.a.
Average weight	800 grams
Max. user weight	100 kg
Ready for use	Yes
Available sizes	3
Machine washable	Yes. Max. 60°C
Drying time	Air dry 5 to 6 hours
Latex free	Yes
Warranty	24 months Lifetime Warranty on SmartCells

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au