



VICAIR ALLROUNDER 02 CUSHION

\$920.00 Original price was: \$920.00.<u>\$828.00</u>Current price is: \$828.00.



SKU: OAL/

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Active as you are, not all of your activities require or allow the use of a wheelchair. In these situations the activity cushion Vicair AllRounder O2 keeps you comfortable and helps to protect your skin. As it is easily buckled onto your body, it follows you around wherever you go.

Feel protected and secure during your sports- and outdoor activities, when you sit in a normal chair or on any hard surface, or when you play on the ground. Afterwards you can simply put the entire AllRounder O2 in the washing machine to clean it up and make it ready for your next adventure together. Experience the freedom!

FEATURES

- Optimal temperature and moisture regulation
- Fully machine washable
- Stays with you during activities outside of the wheelchair
- Skin protection
- Extremely light to support optimal mobility

CHOOSE THE RIGHT SIZE

The activity cushion Vicair AllRounder O2 is available in three sizes: small, medium and large and has an average weight of only 800 grams.

Measure the width of your hips to check which size matches the size of your body. When seated, hold a tape measure horizontally above the widest point of your upper legs, closely against the front of



your hips. Do not follow your leg contours, keep the tape measure straight. Measure from A to B (or B to A). Choose the larger size if the measured size is in between two sizes.

Small < 38 cm / < 15 inch Medium 38 - 42,5 cm / 15- $16\frac{3}{4}$ inch Large 43 - 48 cm / 17 - 19 inch

SPECIFICATIONS

Height n.a.

Average weight 800 grams

Max. user weight 100 kg

Ready for use Yes

Available sizes 3

Machine washable Yes. Max. 60°C

Drying time Air dry 5 to 6 hours

Latex free Yes

Warranty 24 months | Lifetime Warranty on SmartCells