

VICAIR CENTRE RELIEF O2 WHEELCHAIR SEAT CUSHION



SKU: OCR10

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Vicair Centre Relief O2 is the best choice if you need to completely offload the rectum perineum and coccyx area. By simply removing SmartCells from the centre circle compartment you create a non-contact area. This important feature makes a great difference for persons with localised deep pressure ulcers, rectum or tailbone problems or after perineal surgery. The Centre Relief O2 helps you to sit comfortably and safely, even in the most painful situations

The Vicair Centre Relief O2 is a wheelchair cushion that helps prevent pressure ulcers. It is also a high quality addition to a pressure ulcer treatment plan (stage I-IV).

COVER TYPES

Vicair Top Cushion Cover (Qty 2)

- Ideal to help manage one of the main causes of pressure sores
- Standard cover for the Vicair O2 cushions
- Open bottom allows the cushion to breathe continuously
- Optimal ventilation and comfort
- Machine washable

Vicair Comfair Cushion Cover

- Made out of multi-dimensional stretch material
- Ergonomically shaped cover for better envelopment and immersion
- Air transmission and heat dissipation

- Machine washable

Vicair Incotec Cushion Cover

- Usable in all situations where incontinence and hygiene is an issue
- Fluid and incontinence resistant
- Extremely usable for institutional and home use
- Machine washable

SPECIFICATIONS

| | |
|------------------|--|
| Available height | 10 cm |
| Average weight | 750 grams |
| Max. user weight | 250 kg |
| Ready for use | Yes |
| Standard sizes | 4 (Size Chart) |
| Machine washable | Yes. Max. 60°C |
| Drying time | Air dry to 5 to 6 hours |
| Latex free | Yes |
| Warranty | 24 months on inner cover Lifetime Warranty on SmartCells |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au