

VICAIR CENTRE RELIEF O2 WHEELCHAIR SEAT CUSHION

\$933.80 - \$994.86



SKU: OCR10



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Vicair Centre Relief O2 is the best choice if you need to completely offload the rectum perineum and coccyx area. By simply removing SmartCells from the centre circle compartment you create a non-contact area. This important feature makes a great difference for persons with localised deep pressure ulcers, rectum or tailbone problems or after perineal surgery. The Centre Relief O2 helps you to sit comfortably and safely, even in the most painful situations

The Vicair Centre Relief O2 is a wheelchair cushion that helps prevent pressure ulcers. It is also a high quality addition to a pressure ulcer treatment plan (stage I-IV).

COVER TYPES

Vicair Top Cushion Cover (Qty 2)

- Ideal to help manage one of the main causes of pressure sores
- Standard cover for the Vicair O2 cushions
- Open bottom allows the cushion to breathe continuously
- Optimal ventilation and comfort
- Machine washable

Vicair Comfair Cushion Cover

- Made out of multi-dimensional stretch material
- Ergonomically shaped cover for better envelopment and immersion
- Air transmission and heat dissipation

- Machine washable

Vicair Incotec Cushion Cover

- Usable in all situations where incontinence and hygiene is an issue
- Fluid and incontinence resistant
- Extremely usable for institutional and home use
- Machine washable

SPECIFICATIONS

| | |
|------------------|--|
| Available height | 10 cm |
| Average weight | 750 grams |
| Max. user weight | 250 kg |
| Ready for use | Yes |
| Standard sizes | 4 (Size Chart) |
| Machine washable | Yes. Max. 60°C |
| Drying time | Air dry to 5 to 6 hours |
| Latex free | Yes |
| Warranty | 24 months on inner cover Lifetime Warranty on SmartCells |