



# VICAIR CENTRE RELIEF O2 WHEELCHAIR SEAT CUSHION

\$934.00 - \$1,602.00Price range: \$934.00 through \$1,602.00



SKU: OCR10

## **PRODUCT DESCRIPTION**

## **PURPOSE & BENEFITS**

The Vicair Centre Relief O2 is the best choice if you need to completely offload the rectum perineum and coccyx area. By simply removing SmartCells from the centre circle compartment you create a non-contact area. This important feature makes a great difference for persons with localised deep pressure ulcers, rectum or tailbone problems or after perineal surgery. The Centre Relief O2 helps you to sit comfortably and safely, even in the most painful situations

The Vicair Centre Relief O2 is a wheelchair cushion that helps prevent pressure ulcers. It is also a high quality addition to a pressure ulcer treatment plan (stage I-IV).

## **COVER TYPES**

#### Vicair Top Cushion Cover (Qty 2)

- Ideal to help manage one of the main causes of pressure sores
- Standard cover for the Vicair O2 cushions
- Open bottom allows the cushion to breathe continuously
- Optimal ventilation and comfort
- Machine washable

#### **Vicair Comfair Cushion Cover**

- Made out of multi-dimensional stretch material
- Ergonomically shaped cover for better envelopment and immersion
- Air transmission and heat dissipation



• Machine washable

### **Vicair Incotec Cushion Cover**

- Usable in all situations where incontinence and hygiene is an issue
- Fluid and incontinence resistant
- Extremely usable for institutional and home use
- Machine washable

## **SPECIFICATIONS**

Available height 10 cm

Average weight 750 grams

Max. user weight 250 kg

Ready for use Yes

Standard sizes 4 (<u>Size Chart</u>) Machine washable Yes. Max. 60°C

Drying time Air dry to 5 to 6 hours

Latex free Yes

Warranty 24 months on inner cover | Lifetime Warranty on SmartCells